


SE-9409



Introduction

Elastic cord is a very useful item in the physics class. Just a few of the experiments and demonstrations for which you can use it for are described below.

Example experiments

Horizontal waves, mechanically driven

Use a Sine Wave Generator (WA-9867) to drive the cord. As you vary the frequency, the cord will produce excellent wave patterns at different harmonics. Change the tension and examine its effect on the wavelengths.

Waves driven by hand

Secure one end of the cord and set up standing waves. You will want low frequency waves for this, so use a long length of cord and relatively low tension. You can produce standing waves or produce a wave pulse. Give the cord a sharp pulse and watch it travel down the cord and reflect back. Vertical waves may be performed in the same way as horizontal waves.

A very long spring

Use the cord instead of a spring for large demonstrations. The cord will stretch about two times its length. It is almost impossible to over stretch the cord, which is an advantage over springs.

Hang a mass, let it drop, and observe the period of oscillation. Use a Wireless Motion Sensor (PS-3219) to measure displacement as a function of time.

TIP

To make the elastic cord more visible, dye it with fluorescent dye or use a fluorescent pen.

A safe bungee cord jumper

You can simulate bungee jumping in your classroom by attaching a length of the cord to the ceiling and a mass of about 250 g to the end. Keep the length of the bungee cord to about one-third of the height of the room. If the motion can be kept primarily in the vertical plane, you can monitor the results with a Wireless Motion Sensor.

Energy conversion

Hang a mass from the cord, then pull the mass down and release. The amplitude of each oscillation will be less than the previous. By monitoring the decreases in amplitude, the amount of energy lost per oscillation can be calculated. The total energy of a particle in simple harmonic motion is proportional to the square of the amplitude of motion.

Spring constant

Confirm Hooke's law and determine the spring constant of the cord. Measure an unstretched length of a piece of bungee cord. Hang a known mass m from the end of the cord and measure the cord's stretched length. Determine the spring constant using Hooke's law,

$$k = \frac{F}{x}$$

where $F = mg$ and x is the difference between the stretched and unstretched length. Vary the amount of applied mass and measure the corresponding amount of stretch. Plot applied force versus the amount of stretch.

Experiment files


Download one of several student-ready activities from the PASCO Experiment Library. Experiments include editable student handouts and teacher notes. Visit [pasco.com/freelabs/SE-9409](https://www.pasco.com/freelabs/SE-9409).


Specifications and accessories

Visit the product page at [pasco.com/product/SE-9409](https://www.pasco.com/product/SE-9409) to view the specifications and explore accessories. You can also download experiment files and support documents from the product page.

Technical Support

Need more help? Our knowledgeable and friendly Technical Support staff is ready to answer your questions or walk you through any issues.

 **Chat** [pasco.com](https://www.pasco.com)

 **Phone** 1-800-772-8700 x1004 (USA)
+1 916 462 8384 (outside USA)

 **Email** support@pasco.com

Warranty, copyright, and trademarks

Limited warranty

For a description of the product warranty, see the Warranty and Returns page at [pasco.com/legal](https://www.pasco.com/legal).

Copyright

This document is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

Trademarks

PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information, visit [pasco.com/legal](https://www.pasco.com/legal).