PASCO PAScar (Set of 2) ME-6950

Also sold as:

- PAScar Red (ME-6933)
- PAScar Blue (ME-6934)

Each 250-gram polycarbonate cart includes a spring plunger, magnets, and hook and loop pads for collision studies. The PAScars come in red and blue and are compatible with all PASCO Dynamics Tracks and accessories.

Part descriptions



Upper Tie Point
 Use for attaching a string to pull the cart.


2 Accessory Tray

Use to hold Cart Masses (ME-6757A). Use the threaded holes in the tray to attach various cart accessories.

- **3** Hook and Loop Pads
 Use for inelastic collision experiments.
- 4 Plunger Trigger Use to release the spring plunger to launch the cart or for explosion experiments.
- 5 Magnets
 Use for elastic collision experiments
- Use for elastic collision experiments.
- Wheels The wheels are mounted to low-friction ball-breaings. The wheels retract into the body so that they are protected from drops or being stepped on.
- **Or Lower Tie Point**Use for attaching a string to pull the cart.
- 8 Slot
 Use for inserting a Cart Picket Fence (ME-9804).

Plunger operation

The plunger is used to lauch the cart from a stationary object or another PAScar. It can be set at three different positions to launch the cart at three different speeds. The fourth position is used to make the plunger flush with the cart for inelastic collision experiments.



⚠ CAUTION: The plunger comes out rapidly when the trigger is pushed. Do not hold the plunger end of the PAScar next to your eyes to prevent injury.

Launch the cart

Use the plunger to lauch the cart at three different speeds.

- 1. Push the plunger into the PAScar until you hear or feel a "click" (position #1). Repeat to set the plunger at position #2 or #3.
- 2. Place the end of the plunger next to an object.
- 3. Push the plunger trigger to launch the cart.

Set the plunger flush with the cart

The plunger needs to be pushed all the way in (position #4) when using the hook and loop pads for inelastic collision experiments.

- 1. Press and hold down the plunger trigger.
- While holding the trigger, push the plunger all the way into the PAScar until the end of the plunger is flush with the end of the PAScar.
- 3. Release the plunger trigger.

Additional resources

Product information

Visit the product web page at pasco.com/product/me-6950 for additional information including:

- · Specifications
- · Buying Guide
- Experiments
- Documents



Experiments

Print-ready experiment worksheets are available to download from the PASCO website. Go to pasco.com/freelabs and enter ME-6950 in the Part No. field.

Technical Support

Need more help? Our knowledgeable and friendly Technical Support staff is ready to provide assistance with this or any other PASCO product.

Phone (USA) 1-800-772-8700 (Option 4)

Phone (International) +1 916 462 8384

Online pasco.com/support

Warranty, Copyright, and Trademarks

Limited Warranty

For a description of the product warranty, see the Warranty and Returns page at www.pasco.com/legal.

Copyright

This document is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

Trademarks

PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit www.pasco.com/legal.