


# Wireless Optical Dissolved Oxygen Metal Guard

### PS-3604



# Introduction

The metal guard is designed to protect the Wireless Optical Dissolved Oxygen sensor cap. It also enables the probe to sink more rapidly when the probe is immersed in water. A carabiner is included to attach the included 3-meter rope to the sensor.

# To install the guard

- 1. Remove the rubber boot from the sensor probe.
- 2. Put the probe of the sensor into the metal guard.
- 3. Turn the metal guard until snug-tight.




Figure 1. Wireless Optical Dissolved Oxygen sensor with the metal guard installed.

# To attach the rope

- 1. Turn the carabiner nut to open the carabiner.
- 2. Attach the carabiner to the rope.
- 3. Attach the carabiner to the eye bolt on the transparent cover.
- 4. Turn the nut to close the carabiner.




Figure 2. Carabiner attached to the transparent cover eye bolt on the Wireless Optical Dissolved Oxygen sensor.

# **Technical Support**

For assistance with PASCO products, contact PASCO at:

Address:	PASCO scientific 10101 Foothills Blvd. Roseville, CA 95747-7100
Phone:	916-462-8384
Web:	www.pasco.com
Email:	support@pasco.com

### Warranty, Copyright, and Trademarks

**Limited Warranty** For a description of the product warranty, see the PASCO catalog.

**Copyright** This document is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

**Trademarks** PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries. All other brands, products, or service names are or may be trademarks or services of, and are used to identify, products or services of, their respective owners. For more information visit www.pasco.com/legal.