# Instruction Sheet for the PASCO Model OS-8475

# 75 DIODE LASER FOR Iンパー・®OPTICS SYSTEMS

#### Introduction

The OS-8475 Diode Laser for PASCO Optics Systems is designed specifically as an integrated alternate light source for the PASCO OS-8500 Introductory Optics System. This accessory can also be used as an inexpensive light source option for PASCO OS-9250 (advanced) or OS-9251 (intermediate) Optics Systems.

The magnetic steel base of the Diode Laser incorporates an "alignment edge" that, when placed against the alignment rail of your PASCO Optics System, aligns the laser beam to the optical path of your existing system's components. The laser diode assembly can also be rotated 360° adding flexibility to experiment set up. See Figure 1.





Figure 1: OS-8475 Diode Laser Features

## **Equipment**

- Diode Laser for PASCO Optics Systems.
- 9VDC adapter



### Additional equipment recommended:

- PASCO OS-8500 Introductory Optics System
- PASCO OS-9250 Advanced Optics System (without laser)
- PASCO OS-9251 Intermediate Optics System (without laser)

## **Replacement parts:**

Description PASCO part no.

thumbscrew assembly 617-008 Optics System mount 648-06146

♦ CAUTION - use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

LASER LIGHT - DO NOT STARE INTO BEAM.

#### **CLASS 2 LASER PRODUCT**

This Laser Product is designated as Class 2 during all procedures of operation.

Do not point the laser or allow the laser to be directed or reflected toward people or animals.





## **Specifications**

Classification	Class II Laser
Wavelength	532 nm (green) 650 nm (red)
Maximum Output	<1 mW
Divergence	<2 mrd
Beam Diameter at aparture	3 mm

### **Assembly**

- 1. Determine the configuration of your experiment and the location in which the Diode Laser will be set up and adjust the orientation of the Diode Laser accordingly. Make sure the face of the laser diode assembly is relatively perpendicular to the surface of your optics bench. This will minimize the amount of fine adjustment you may need to perform in step 4. Tighten the thumbscrew to secure the laser diode module.
- 2. Place the "alignment edge" of the Diode Laser against the alignment rail of your PASCO optics bench. See Figures 2 and 3.
- 3. Plug the 9VDC adapter phone plug into the back of the laser and turn on the laser.
- 4. If necessary, use the adjustment screws to align the laser beam from left-to-right and up-and down. See Figure 4.
- 5. When using with the Laser Slits, OS-8529, slide the Laser Diode to the top of the slots on the mount




Figure 2: Laser Diode Setup foPASCO Introductory
Optics bench




Figure 3: Laser Diode Setup for PASCO Advanced Optics bench




Figure 4: Laser Beam Fine Adjustment

2 12 H 40 ®

## **Technical Support**

For assistance with any PASCO product, contact PASCO at:

Address: PASCO scientific

10101 Foothills Blvd. Roseville, CA 95747-7100

Phone: 916-786-3800 (worldwide)

800-772-8700 (U.S.)

Fax: (916) 786-7565 Web: www.pasco.com Email: support@pasco.com

#### **Limited Warranty**

For a description of the product warranty, see the PASCO catalog.

NOTE: NO service or maintenance is allowed, on this product, by the customer. Return unit to the factory for service or repair. Unit is NOT to be opened or modified by the Customer.

#### Copyright

The PASCO scientific 012-09867A *Basic Optics Diode Laser Instruction Sheet* is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

#### **Trademarks**

PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit www.pasco.com/legal.



# Notes:

4 II:E∀€•®