



Dynamics Track End Stop (ME-8971)

Introduction



The Dynamics Track End Stops are intended to be used with PASCO tracks and carts. The two included end stops are compatible with any aluminum or plastic PASCO track, and with any PASCO cart. One side of each end stop contains magnets to repel the magnetic bumper of a cart. The stud on top of each end stop can be used for attaching a string or a spring to the end stop. The hole in the end stop is meant to be used with a Spring Cart Launcher (ME-6843).

Features



- 1 T-slot tab
- 2 Thumbscrew
- 3 Magnets
- 4 String and spring attachment stud
- 5 Hole for Spring Cart Launcher

Installation

Attach each end stop to the side of the tracks using these steps:




NOTE: If you are using an older PASCO track, see the **Using the End Stops with older tracks** section before installing the End Stops.

1. Loosen the square nut on the end of the thumbscrew until it is able to fit into the T-slot on the side of the track.
2. Slide the tab and square nut into the T-slot.
3. Position the end stop at the desired distance from the end of the track.
4. Tighten the thumbscrew to secure the end stop in place.

Using the End Stops with older tracks

Older PASCO tracks feature T-slots with less clearance than the T-slots of current models. If you have an older PASCO track, install the included washers on the End Stop thumbscrews, as shown below, to prevent the thumbscrews from "bottoming out".



Experiment files

Download one of several student-ready activities from the PASCO Experiment Library. Experiments include editable student handouts and teacher notes. Visit pasco.com/freelabs/ME-8971.

Technical support

Need more help? Our knowledgeable and friendly Technical Support staff is ready to answer your questions or walk you through any issues.

- ☐ Chat pasco.com
- Phone 1-800-772-8700 x1004 (USA)
+1 916 462 8384 (outside USA)
- ☐ Email support@pasco.com

Regulatory information

Limited warranty

For a description of the product warranty, see the Warranty and Returns page at www.pasco.com/legal.

Copyright

This document is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

Trademarks

PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit www.pasco.com/legal.