

Flex Rod

ME-8978A



Included Equipment 1. Flexible Tubing 2. Rod Clamp 3. Lock Nut (2) 4. Three-finger Clamp 5. Base Adapter Rod 6. Double Clamp

Introduction

The PASCO ME-8978A Flex Rod is a unique variation on the traditional lab stand. Use it with a base or table clamp (sold separately) to support sensors, photogates, pulleys, diode lasers, test tubes, and other lightweight equipment. The flexible "goose-neck" tubing allows easily and accurate positioning by bending the rod.

Thumbscrew (2)

Sensors and other equipment can be screwed directly to Flex Rod's 1/4-20 threaded studs or held in the included rod clamp, double clamp, or three-finger clamp. The rod clamp and double clamp can also be screwed directly to Flex Rod's threaded studs. The rod clamp and double clamp hold cylindrical items up to 12 millimeters (mm) in diameter. The three finger clamp spans 62 mm when fully open.

Compatible Bases

The assembled Flex Rod can be mounted on any base or table clamp that accepts an unthreaded 12.7 mm (0.5 inch) diameter rod or a 1/2-13 threaded rod.

The table below lists compatible bases and table clamps available from PASCO.

Compatible Base	Part Number
Large Rod Base	ME-8735
Small "A" Base	ME-8976
Large Base and Support Rod	ME-9355
Large Table Clamp	ME-9472
Universal Table Clamp	ME-9376B

800-772-8700 www.pasco.com

Model No. ME-8978A Set-up

Flexible Tubing

Lock Nut

Rod


Base Adapter


Clamp or

screw into


Set-up

- Screw one end of the flexible tubing through a lock nut into the mating hole of the Base Adapter Rod.
- 2. Secure the Base Adapter Rod in a base or table clamp either by clamping it or screwing it into a 1/2-13 threaded hole.
- 3. Screw a second lock nut onto the threaded stud on the free end of the Flexible Tubing.






5. To attach the included rod clamp or double clamp, screw it loosely onto the threaded stud; then tighten the lock nut against it.



6. To attach the included three-finger clamp, first attach a rod clamp or double clamp to the flexible tubing. Mount the three-finger clamp in the rod clamp and tighten the thumbscrew to hold it in place.



Technical Support

For assistance with any PASCO product, contact PASCO at:

Address: PASCO scientific

10101 Foothills Blvd. Roseville, CA 95747-7100

Phone: 916-462-8384 (worldwide)

877-373-0300 (U.S.)

Web: www.pasco.com
Email: support@pasco.com

Limited Warranty

For a description of the product warranty, see the PASCO catalog.

Copyright

The PASCO scientific 012-10009B Flex Rod Instruction Sheet is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

Trademarks

PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit www.pasco.com/legal.