

Complete Microwave Optics System

WA-9316A

Experiment Guide



Important Note

• CAUTION — Under some circumstances, microwaves can interfere with electronic medical devices. If you use a pacemaker, or other electronic medical device, check with your doctor or the manufacturer to be certain that low power microwaves at a frequency of 10.525 GHz will not interfere with its operation.



Table Of Contents

Introduction 3
Microwave Optics System Equipment 3
Technical Support 4
Experiments
Experiment 1: Introduction to the Microwave Optics System 5
Experiment 2: Reflection 9
Experiment 3: Standing Waves - Measuring Wavelengths11
Experiment 4: Refraction Through a Prism15
Experiment 5: Polarization17
Experiment 6: Double-Slit Interference19
Experiment 7: Lloyds Mirror21
Experiment 8: Fabry-Perot Interferometer23
Experiment 9: Michelson Interferometer25
Experiment 10: Fiber Optics27
Experiment 11: Brewster's Angle28
Experiment 12: Bragg Diffraction31
(Note that Experiment 3 can be done with the WA-9319A Microwave Detector Probe and that Experiment 11 and Experiment 12 use equipment from the WA-9315 Microwave Accessory Package.)
Teacher's Guide

Introduction

The experiments are designed to give a thorough introduction to wave theory using the Microwave Optics System. Basic principles are only briefly discussed in each experiment because it is assumed that the students have the appropriate theoretical background. These experiments range from quantitative investigation of reflection and refraction to microwave models of the Michelson and Fabry-Perot interferometers. Experiments for investigating Bragg diffraction and Brewster's angle are included, but require the Microwave Accessory Package.

The experiments are written in a worksheet format. Fee free to copy the experiment write-ups for use in your lab.

Microwave Optics System Equipment

The Complete Microwave Optics System consists of the Microwave Optics System (WA-9314C) and the Microwave Accessory Package (WA-9315). The Microwave Optics System comes with a 2.85 centimeter wavelength microwave transmitter and a receiver with variable amplification from 1X to 30X. All the accessory equipment needed to investigate a variety of wave phenomena is also included.

Microwave Optics System (WA-9314C)	Microwave Optics System (WA-9314C)
Microwave Receiver WA-9800	Metal Reflectors (2)
Microwave Transmitter WA-9801	Partial Reflectors (2)
Mounting Stands (2) WA-9802	Polarizers (2)
Power Adapters (2)	Slit Extender Arm
Goniometer	Narrow Slit Spacer
Fixed Arm Assembly	Wide Slit Spacer
Rotating Table	Ethafoam Prism Mold
Component Holders (2)	Styrene Pellets (1 jar)
Rotating Component Holder	Tubular Plastic Bags (4)

Microwave Accessory Package (WA-9315)
Cubic Lattice* (Bragg's crystal model)
Polyethylene Panel

*The Cubic Lattice contains 100 metal spheres in a 5 by 5 by 4 array.

Assembling Equipment for Experiments

The arms of the Goniometer or the Fixed Arm Assembly slide through the holes in the Component Holders as shown. Make sure the magnetic strip on the bottom of the arm grips the base of the carriage. To adjust the position of the holders, just slide them along the arms. Attach the mounting stands of the microwave Transmitter and Receiver to the arms of the Goniometer or Fixed Arm Assembly in the same manner.

For most experiments it is advantageous to attach the Transmitter to the long arm of the Goniometer and the Receiver to the shorter, rotatable arm. This maintains a fixed relationship between the microwave beam and components mounted on the long arm (or on the degree plate) of the Goniometer. In turn the Receiver moves easily to sample the output.



Reflectors, Partial Reflectors, Polarizers, Slit Spacers, and the Slit Extender Arm all attach magnetically to the Component Holders. The metric scale along the Goniometer and Fixed Arm Assembly arms and the degree plate at the junction of the Goniometer arms allow easy measurement of component placement.

- It is recommended that the apparatus be mounted on a CLEAN, SMOOTH table.
- Before setting up the equipment, brush off any material—particularly metal chips—that might have adhered to the magnetic strips and magnetic pads.

Maintenance

Regular maintenance for this equipment is minimal. When not in use, the equipment can be stored in the shipping box. Inspect the magnetic strips and magnetic pads and remove anything - especially bits of metal - that adhere to them.

If problems arise with the equipment DO NOT attempt to fix this equipment yourself. Contact PASCO scientific. (See the Technical Support information.)

Technical Support

For assistance with any PASCO product, contact PASCO at:

Address: PASCO scientific 10101 Foothills Blvd. Roseville, CA 95747-7100

Phone: +1 916 786 3800 (worldwide) 800-772-8700 (U.S.) E-mail: support@pasco.com

Web www.pasco.com

For the latest information about the Complete Microwave Optics System or its accessories,

go to the PASCO web site at www.pasco.com and enter the model number or the product name in the search window.

Limited Warranty For a description of the product warranty, see the PASCO catalog. Copyright The PASCO scientific *Experiment Guide* is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited. Trademarks PASCO, PASCO Capstone, PASPORT, SPARK Science Learning System, SPARK SLS, and SPARKvue are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries. For more information visit www.pasco.com/legal.

Product End of Life Disposal Instructions:

This electronic product is subject to disposal and recycling regulations that vary by country and region. It is your responsibility to recycle your electronic equipment per your local environmental laws and regulations to ensure that it will be recycled in a manner that protects human health and the environment. To find out where you can drop off your waste equipment for recycling, please contact your local waste recycle/disposal service, or the place where you purchased the product.

The European Union WEEE (Waste Electronic and Electrical Equipment) symbol (to the right) and on the product or its packaging indicates that this product must not be disposed of in a standard waste container.



Experiment 1: Introduction to the Microwave Optics System

Equipment Needed:

Item	Item	
Transmitter	Goniometer	
Receiver	Reflector	

Purpose

This activity gives an introduction to the Microwave Optics System which may prove helpful in learning to use the equipment effectively and in understanding the significance of measurements made with this equipment.

Procedure

- 1. Arrange the Transmitter and Receiver on the Goniometer as shown in Figure 1 with the Transmitter attached to the fixed arm. Be sure to adjust both Transmitter and Receiver to the same polarity — the horns should have the same orientation, as shown.
- 2. Plug in the Transmitter and the Receiver. Turn the INTENSITY selection switch on the Receiver from OFF to 10X. (The LEDs should light up on both units.)
- 3. Adjust the Transmitter and Receiver so the distance between the source diode in the Transmitter and the detector diode in the Receiver (the distance labeled **R** in Figure 1) is 40 cm (see Figure 1.2 for location of points of transmission and reception). The diodes are at the locations directly above the points marked "T" and "R" on the stands. Adjust the INTENSITY and VARIABLE SENSITIVITY knobs on the Receiver so that the meter reads 1.0 (full scale).



4. Set the distance **R** to each of the values shown in Table 1. Record the meter reading for each value of **R**. (Do not adjust the Receiver controls between measurements.) After making the measurements, perform the calculations shown in the table.

Table 1.1: Data	а
-----------------	---

R (cm)	Meter Reading (M)	M X R (cm)	M X R ² (cm ²)
40			
50			
60			
70			
80			
90			
100			

Experiment Guide

- 5. Set **R** to some value between 70 and 90 cm. While watching the meter, slowly decrease the distance between the Transmitter and Receiver. Does the meter deflection increase steadily as the distance decreases?
- 6. Set R to between 50 and 90 cm. Move a Reflector, its plane parallel to the axis of the microwave beam, toward and away from the beam axis, as shown in Figure 1.3. Observe the meter readings. Can you explain your observations in steps 5 and 6? (Don't worry if you can't; you will



have a chance to investigate these phenomena more closely in other experiments.) For now just be aware of the following:

IMPORTANT: Reflections from nearby objects, including the table top, can affect the results of your microwave experiments. To reduce the effects of extraneous reflections, keep your experiment table clear of all objects, especially metal objects, other than those components required for the current experiment.

7. Loosen the hand screw on the back of the Receiver and rotate the Receiver as shown in Figure 1.4. This varies the polarity of maximum detection. (Look into the receiver horn and notice the alignment of the detector diode.) Observe the meter readings through a full 360 degree rotation of the horn. A small mirror may be helpful to view the meter reading as the receiver is turned. At what polarity does the Receiver detect no signal?

NOTE: Try rotating the Transmitter horn as well. When finished, reset the Transmitter and Receiver so their polarities match (e.g., both horns are horizontal or both horns are vertical).

8. Position the Transmitter so the output surface of the horn is centered directly over the center of the Degree Plate of the Goniometer arm (see Figure 1.5). With the Receiver directly facing the Transmitter and as far back on the Goniometer arm as possible, adjust the Receiver controls for a meter reading of 1.0. Then rotate the rotatable arm of the Goniometer as shown in the figure. Set the angle of rotation (measured relative to the 180-degree point on the degree scale) to each of the values shown in Table 2, and record the meter reading at each setting



Table 1.2: Data

Angle of Receiver	Meter Reading	Angle of Receiver	Meter Reading	Angle of Receiver	Meter Reading
90°		160		230	
100		170		240	
110		180		250	
120		190		260	
130		200		270	
140		210			
150		220			



Questions

- 1. The electric field of an electromagnetic wave is inversely proportional to the distance from the wave source, (i.e., E = 1/R). Use your data from step 4 of the experiment to determine if the meter reading of the Receiver is directly proportional to the electric field of the wave.
- 2. The intensity of an electromagnetic wave is inversely proportional to the square of the distance from the wave source (i.e., $I = 1/R^2$). Use your data from step 4 of the experiment to determine if the meter reading of the Receiver is directly proportional to the intensity of the wave.
- **3.** Considering your results in step 7, to what extent can the Transmitter output be considered a spherical wave? A plane wave?

Experiment 2: Reflection

Equipment Needed:

Item	Item	
Transmitter	Goniometer	
Receiver	Metal Reflector	
Rotating Component Holder		

Purpose

This experiment explores the relationship of reflected microwaves to incident microwaves when the incident microwaves encounter a flat reflector at some angle measured from the normal of the reflector.

Procedure

- 1. Arrange the equipment as shown in Figure 2.1 with the Transmitter attached to the fixed arm of the Goniometer. Be sure to adjust the Transmitter and Receiver to the same polarity; the horns should have the same orientation as shown.
- **2.** Plug in the Transmitter and turn the Receiver INTENSITY selection switch to 30X.
- **3.** The angle between the incident wave from the Transmitter and a line normal to the plane of the Reflector is called the Angle of Incidence (see Figure 2.2). Adjust the Rotating Component Holder so that the Angle of Incidence equals 45°.
- 4. Without moving the Transmitter or the Reflector, rotate the movable arm of the Goniometer until the meter reading is a maximum. The angle between the axis of the Receiver horn and a line normal to the plane of the Reflector is called the Angle of Reflection.



5. Measure and record the angle of reflection for each of the angles of incidence shown in Table 2.1.

Table 2.1: Data			
Angle of Incidence	Angle of Reflection		
20°			
30°			
40°			
50°			
60°			
70°			
80°			
90°			

NOTE: At various angle settings the Receiver will detect both the reflected wave and the wave coming directly from the Transmitter, thus giving misleading results. Determine the angles for which this is true and mark the data collected at these angles with an asterisk ("*").

Questions

- 1. What relationship holds between the angle of incidence and the angle of reflection? Does this relationship hold for all angles of incidence?
- 2. In measuring the angle of reflection, you measured the angle at which a maximum meter reading was found. Can you explain why some of the wave reflected into different angles? How does this affect your answer to question 1?
- **3.** Ideally you would perform this experiment with a perfect plane wave, so that all the Transmitter radiation strikes the Reflector at the same angle of incidence. Is the microwave from the Transmitter a perfect plane wave (see Experiment 1, step 7)? Would you expect different results if it were a perfect plane wave? Explain.

Questions for Additional Experimentation

- 1. How does reflection affect the intensity of the microwave? Is all the energy of the wave striking the Reflector reflected? Does the intensity of the reflected signal vary with the angle of incidence?
- 2. Metal is a good reflector of microwaves. Investigate the reflective properties of other materials. How well do they reflect? Does some of the energy pass through the material? Does the material absorb some of it? Compare the reflective properties of conductive and non-conductive materials.

Experiment 3: Standing Waves – Measuring Wavelength

Equipment Needed:

Item	Item
Transmitter	Goniometer
Receiver	Reflector
Component Holder (2)	Optional: Microwave Detector Probe

Introduction

When two electromagnetic waves meet in space, they superpose. The superposition principle describes that the total electric field at any point is the sum of the electric fields created by both waves at that point. If the two waves travel at the same frequency but in opposite direction, they form a standing wave. Nodes appear where the fields of the two waves cancel and antinodes appear where the superposed fields oscillate between a maximum and a minimum. The distance between nodes in the standing wave pattern is just one-half (1/2) the wavelength, λ , of the two waves.

Purpose

In this experiment, you will reflect the wave from the Transmitter back upon itself, creating a standing wave pattern. By measuring the distance between nodes in the pattern and multiplying by two, you can determine the wavelength, λ , of the microwave radiation.

Procedure

Method A

1. Set up the equipment as shown in Figure 3.1. Adjust the Receiver controls to get a full-scale meter reading with the Transmitter and Receiver as close together as possible. Slowly move the Receiver along the Goniometer arm away from the Transmitter. How does this motion effect the meter reading?

NOTE: The microwave horns are not perfect collectors of microwave radiation. Instead, they act as partial reflectors, so that the radiation from the Transmitter reflects back and forth between the Transmitter and Reflector horns, diminishing in amplitude at each pass. However, if the distance



between the Transmitter and Receiver diodes is equal to $n\lambda/2$, (where **n** is an integer and λ is the wavelength of the radiation), then all the multiply-reflected waves entering the Receiver horn will be in phase with the primary transmitted wave. When this occurs, the meter reading will be a maximum. (The distance between adjacent positions in order to see a maximum is therefore $\lambda/2$.)

- 2. Slide the Receiver one or two centimeters along the Goniometer arm to obtain a maximum meter reading. Record the Receiver position along the metric scale of the Goniometer arm.
- Initial Position of Receiver = ______
- 3. While watching the meter, slide the Receiver away from the Transmitter. Do not stop until the Receiver passed through at least 10 positions at which you see a minimum meter reading and it returned to a position where the reading is a maximum. Record the new position of the Receiver and the number of minima that were traversed.
- Minima Traversed=_____.
- Final Receiver Position = _____.
- 4. Use the data you have collected to calculate the wavelength of the microwave radiation.

Experiment Guide

- λ = _____.
- 5. Repeat your measurements and recalculate l.
- Initial Position of Receiver = _____.
- Minima Traversed = _____.
- Final Receiver Position = _____.
- λ = _____

Question

Use the relationship velocity = λv to calculate the frequency, v, of the microwave signal (assuming velocity of propagation in air is 3 x 10⁸ m/s).

(v = the expected frequency of the microwave radiation - 10.525 GHz).

NOTE: This experiment can also be performed using the PASCO Microwave Detector Probe (Model ME-9319A), as described in Method B below. Method A can be used for those without a probe. However, in Method A, the wavelength, λ , can not be measured *directly* from the standing wave pattern. In Method B, the wavelength can be measured directly from the standing wave pattern.

Method B

- 1. Arrange the equipment as shown in Figure 3.1.
- 2. Plug the Detector Probe into the EXTERNAL ANTENNA connector on the Receiver. Face the Receiver horn directly *away* from the Transmitter so that none of the microwave signal enters the horn. Adjust the Receiver controls as needed to get a strong meter reading.
- 3. Slide the Probe along the Goniometer arm (no more than a centimeter or two) until the meter shows a maximum read-



ing. Then slide the Reflector (again, no more than a centimeter or two) to find a maximum meter reading. Continue making slight adjustments to the Probe and Reflector positions until the meter reading is as high as possible.

- 4. Now find a node of the standing wave pattern by adjusting the Probe until the meter reading is a minimum. Record the Probe Position along the metric scale on the Goniometer arm.
- Initial Probe Position =
- 5. While watching the meter, slide the Probe along the Goniometer arm until the Probe has passed through at least 10 antinodes and returned to a node. Record the new position of the Probe and the number of antinodes that were traversed.
- Antinodes Traversed=_____.
- Final Probe Position = _____.
- 6. Use your data to calculate λ , the wavelength of the microwave radiation.

λ=____.

- 7. Repeat your measurements and recalculate λ .
- Initial Probe Position =_____.
- Antinodes Traversed =_____.

- Final Probe Position =_____.
- λ =_____.

Question

Use the relationship velocity = λv to calculate the frequency, v, of the microwave signal (assuming velocity of propagation in air is 3 x 10⁸ m/s).

(v = the expected frequency of the microwave radiation - 10.525 GHz).

Experiment 4: Refraction Through a Prism

Equipment Needed:

Item	Item
Transmitter	Goniometer
Receiver	Rotating Table
Prism Mold with Styrene Pellets	Protractor

Introduction

An electromagnetic wave usually travels in a straight line. As it crosses a boundary between two different media, however, the direction of propagation of the wave changes. This change in direction is called Refraction, and it is summarized by a mathematical relationship known as the **Law of Refraction** (otherwise known as **Snell's Law**):

$$\mathbf{n}_1 \mathbf{sin} \theta_1 = \mathbf{n}_2 \mathbf{sin} \theta_2$$

where θ_1 is the angle between the direction of propagation of the incident wave and the normal to the boundary between the two media, and θ_2 is the corresponding angle for the refracted wave (see Figure 4.1). Every material can be described by a number **n**, called its **Index of Refraction**. This number indicates the ratio between the speed of electromagnetic waves in vacuum and the speed of electromagnetic waves in the material, also called the medium. In general, the media on either side of a boundary will have different indices of refraction. Here they are labeled **n**₁ and **n**₂. It is the dif-



ference between indices of refraction (and the difference between wave velocities that this implies) which causes "bending", or refraction of a wave as it crosses the boundary between two distinct media.

Purpose

In this experiment, you will use the law of refraction to measure the index of refraction for styrene pellets.

Procedure

- 1. Arrange the equipment as shown in Figure 4.2. Rotate the empty prism mold and see how it effects the incident wave. Does it reflect, refract, or absorb the wave?
- 2. Fill the prism mold with the styrene pellets. To simplify the calculations, align the face of the prism that is nearest to the Transmitter so it is perpendicular to the incident microwave beam.
- 3. Rotate the movable arm of the Goniometer and locate the angle θ at which the refracted signal is a maximum.



NOTE: θ is just the angle that you read directly from the Degree Scale of the Goniometer.

• θ = _____.

Experiment Guide

- 4. Using the diagram shown in Figure 4.3, determine θ_1 and use your value of θ to determine θ_2 . (You will need to use a protractor to measure the Prism angles.)
- $\theta_1 =$ _____.
- $\theta_2 =$ _____.
- 5. Plug these values into the Law of Refraction to determine the value of n_1/n_2 .
- $\mathbf{n}_1/\mathbf{n}_2 =$ _____.



Questions

- 1. In the diagram of Figure 4.3, the assumption is made that the wave is unrefracted when it strikes the first side of the prism (at an angle of incidence of 0°). Is this a valid assumption?
- 2. Using this apparatus, how might you verify that the index of refraction for air is equal to 1.00.
- 3. Would you expect the refraction index of the styrene pellets in the prism mold to be the same as for a solid styrene prism?



Experiment 5: Polarization

Equipment Needed:

Item	Item	
Transmitter	Goniometer	
Receiver	Polarizer (1)	
Component Holder (2)		



Introduction

The microwave radiation from the Transmitter is linearly polarized along the Transmitter diode axis (i.e., as the radiation propagates through space, its electric field remains aligned with the axis of the diode). If the Transmitter diode were aligned vertically, the electric field of the transmitted wave would be vertically polarized, as shown in Figure 5.1. If the detector diode were at an angle θ to the Transmitter diode, as shown in Figure 5.2, it would only detect the component of the incident electric field that was aligned along its axis.

Purpose

In this experiment you will investigate the phenomenon of polarization and discover how a polarizer can be used to alter the polarization of microwave radiation.

Procedure

- 1. Arrange the equipment as shown in Figure 5.3 and adjust the Receiver controls for nearly full-scale meter deflection.
- 2. Loosen the hand screw on the back of the Receiver and rotate the Receiver in increments of ten degrees. At each rotational position, record the meter reading in Table 5.1.
- 3. What happens to the meter readings if you continue to rotate the Receiver beyond 180-degrees



Table 5.1: Data

Angle of Receiver	Meter Reading	Angle of Receiver	Meter Reading	Angle of Receiver	Meter Reading
0°		70°		140°	
10°		80°		150°	
20°		90°		160°	
30°		100°		170°	
40°		110°		180°	
50°		120°			
60°		130°			

- 4. Set up the equipment as shown in Figure 5.4. Reset the Receivers angle to 0-degrees (the horns should be oriented as shown with the longer side horizontal).
- 5. Record the meter reading when the Polarizer is aligned at 0, 22.5, 45, 67.5 and 90-degrees with respect to the horizontal.
- 6. Remove the Polarizer slits. Rotate the Receiver so the axis of its horn is at right angles to that of the Transmitter. Record the meter reading. Then replace the Polarizer slits and record the meter readings with the Polarizer slits horizontal, vertical, and at 45-degrees.

Angle of Polarizer	Meter Reading
0° (horizontal)	
22.5°	
45°	
67.5°	
90° (vertical)	



Angle of Slits	Meter Reading
Horizontal	
Vertical	
45°	

Questions

- 1. If the Receiver meter reading (M) were directly proportional to the electric field component (E) along its axis, the meter would read the relationship $\mathbf{M} = \mathbf{M}_0 \cos\theta$ (where θ is the angle between the detector and Transmitter diodes and \mathbf{M}_0 is the meter reading when $\theta = 0$). (See Figure 5.2). Graph your data from step 2 of the experiment. On the same graph, plot the relationship $\mathbf{M}_0 \cos\theta$. Compare the two graphs.
- 2. The intensity of a linearly polarized electromagnetic wave is directly proportional to the square of the electric field (e.g., $I = kE^2$). If the Receiver's meter reading was directly proportional to the incident microwave's intensity, the meter would read the relationship $M = M_0 \cos^2 \theta$. Plot this relationship on your graph from question 1. Based on your graphs, discuss the relationship between the meter reading of the Receiver and the polarization and magnitude of the incident microwave.
- 3. Based on your data from step 5, how does the Polarizer affect the incident microwave?
- 4. Can you explain the results of step 6 of the experiment. How can the insertion of an additional polarizer increase the signal level at the detector? (HINT: Construct a diagram like that shown in Figure 5.2 showing (1) the wave from the Transmitter; (2) the wave after it passes through the Polarizer; and (3) the component detected at the detector diode.)

Experiment 6: Double-Slit Interference

Equipment Needed:

Item	Item
Transmitter	Goniometer
Receiver	Rotating Component Holder)
Slit Extender Arm	Metal Reflectors (2)
Wide Slit Spacer	Narrow Slit Spacer

Introduction

Experiment 3 showed that two waves moving in opposite directions can superpose to create a standing wave pattern. A somewhat similar phenomenon occurs when an electromagnetic wave passes through a two-slit aperture. The wave *diffracts* into two waves which superpose in the space beyond the apertures. Similar to the standing wave pattern, there are points in space where maxima are formed and others where minima are formed.

With a double slit aperture, the intensity of the wave beyond the aperture will vary depending on the angle of detection. For two thin slits separated by a distance d, maxima will be found at angles such that



d $\sin\theta = n\lambda$. (where θ = the angle of detection, λ = the wavelength of the incident radiation, and **n** is any integer; see Figure 6.1). Refer to a textbook for more information about the nature of the double-slit diffraction pattern.

Purpose

In this experiment you will investigate the phenomena of diffraction and interference.

Procedure

- 1. Arrange the equipment as shown in Figure 6.2. Use the Slit Extender Arm, two Reflectors, and the Narrow Slit Spacer to construct the double slit. (A slit width of about 1.5 cm is recommended.) Be precise with the alignment of the slit and make the setup as symmetrical as possible.
- 2. Adjust the Transmitter and Receiver for vertical polarization (0°) and adjust the Receiver controls to give a full-scale reading at the lowest possible amplification.
- **3.** Rotate the rotatable Goniometer arm (on which the Receiver rests) slowly about its axis. Observe the meter readings.



- 4. Reset the Goniometer arm so the Receiver directly faces the Transmitter. Adjust the Receiver controls to obtain a meter reading of 1.0. Now set the angle θ to each of the values shown in Table 6.1. At each setting record the meter reading in the table. (In places where the meter reading changes significantly between angle settings, you may find it useful to investigate the signal level at intermediate angles.)
- 5. Keep the slit widths the same, but change the distance between the slits by using the Wide Slit Spacer instead of the Narrow Slit Spacer. Because the Wide Slit Space is 50% wider than the Narrow Slit Spacer (9 cm vs 6 cm) move the Trans-

mitter back 50% so that the microwave radiation at the slits will have the same relative intensity. Repeat the measurements. (You may want to try other slit spacings as well.)

Table 6.1: Data					
Angle	Meter Reading	Angle	Meter Reading	Angle	Meter Reading
0°		30°		60°	
5°		35°		65°	
10°		40°		70°	
15°		45°		75°	
20°		50°		80°	
25°		55°		86°	

Questions

- 1. From your data, plot a graph of meter reading versus q. Identify the angles at which the maxima and minima of the interference pattern occur.
- 2. Calculate the angles at which you would expect the maxima and minima to occur in a standard two-slit diffraction pattern maxima occur wherever $\mathbf{d} \sin \theta = \mathbf{n} \lambda$, and minima occur wherever $\mathbf{d} \sin \theta = \mathbf{n} \lambda/2$. (Check your textbook for the derivation of these equations, and use the wavelength measured in experiment 3.) How does this compare with the locations of your observed maxima and minima? Can you explain any discrepancies? (What assumptions are made in the derivations of the formulas and to what extent are they met in this experiment?)
- **3.** Can you explain the relative drop in intensity for higher order maxima? Consider the single-slit diffraction pattern created by each slit. How do these single slit patterns affect the overall interference pattern?

NOTE:

- Wavelength at 10.525 GHz = 2.85 cm.
- The experimenter's body position may affect the results.

Experiment 7: Lloyd's Mirror

Equipment Needed:

Item	Item
Transmitter	Goniometer
Receiver	Component Holder)
Fixed Arm Assembly	Reflector (1)
Meter Stick	

Introduction

In earlier experiments, such as 3 and 6, you observed how a single electromagnetic wave can be diffracted into two waves and, when the two components join back together, they form an interference pattern. Lloyd's Mirror is another example of this phenomenon. Just as with the other interference patterns you have seen, this interference pattern provides a convenient method for measuring the wavelength of the radiation.

Figure 7.1 is a diagram for Lloyd's Mirror. An electromagnetic wave from point source **A** is detected at point **C**. Some of the electromagnetic wave, of course, propagates directly between point **A** and **C**, but some reaches C after being reflected at point **B**. A maximum signal



will be detected when the two waves reach the detector in phase. Assuming that the diagram shows a setup for a maximum signal, another maximum will be found when the Reflector is moved back so the path length of the reflected beam is $AB + BC + \lambda$. where λ is the wavelength of the electromagnetic wave,

Purpose

In this experiment you will investigate the phenomenon of an interference pattern formed by a Lloyd's Mirror.

Procedure

- Arrange the equipment as shown in Figure 7.2. For best results, the Transmitter and Receiver should be as far apart as possible. Be sure the Receiver and Transmitter are equidistant (d₁) from the center of the Goniometer degree plate and that the horns are directly facing each other. (See Figure 7.3 for location of effective points of transmission and reception). Also be sure that the surface of the Reflector is parallel to the axis of the Transmitter and Receiver horns.
- 2. While watching the meter on the Receiver, slowly slide the Reflector away from the Degree Plate. Notice how the meter reading passes through a series of minima and maxima.
- **3.** Find the Reflector position closest to the degree plate which produces a minimum meter reading.



4. Measure and record \mathbf{h}_1 , the distance between the center of the degree plate and the surface of the Reflector.

^{5.} Slowly slide the Reflector away from the degree plate until the meter reading passes through a maximum and returns to a new minimum. Measure and record \mathbf{h}_2 , the new distance between the center of the degree plate and the surface of the Reflector.



[•] **h**₁ = _____

Experiment Guide

- **h**₂ = _____.
- 6. Measure d_1 the distance between the center of the degree scale and the Transmitter diode.
- **d**₁ = _____.
- 7. Use your collected data to calculate λ , the wavelength of the microwave radiation.
- λ = _____.
- 8. Change the distance between the Transmitter and Receiver and repeat your measurements.
- **h**₁ = _____.
- **h**₂ = _____.
- d₁ = _____.
- λ = _____.



Question

- What is the advantage in having the effective transmission and reception points equidistant from the center of the degree plate in this experiment?
- **NOTE**: Don't stand in front of the apparatus while conducting the experiment. Your body acts as a reflector. Therefore, try to stand to one side behind the plane of the receiver horn.



Experiment 8: Fabry-Perot Interferometer

Equipment Needed:

Item	ltem
Transmitter	Goniometer
Receiver	Component Holders (2)
Partial Reflectors (2)	

Introduction

When an electromagnetic wave encounters a partial reflector, part of the wave reflects and part of the wave transmits through the partial reflector. A Fabry-Perot Interferometer consists of two parallel partial reflectors positioned between a wave source and a detector (see Figure 8.1).

The wave from the source reflects back and forth between the two partial reflectors. However, with each pass, some of the radiation passes through to the detector. If the distance between the partial reflectors is equal to $n\lambda/2$, where λ is the wavelength of the radiation and **n** is an integer, then all the waves passing through to the detector at any instant will be in phase. In this case, a maximum signal will be detected by the Receiver. If the distance between the partial reflectors is not a multiple of $\lambda/2$, then some degree of destructive interference will occur, and the signal will not be a maximum.

Purpose

In this experiment you will investigate the Fabry-Perot Interferometer.

Procedure

- 1. Arrange the equipment as shown in Figure 8.1. Plug in the equipment and adjust the Receiver controls for an easily readable signal.
- 2. Adjust the distance between the Partial Reflectors and observe the relative minima and maxima.
- 3. Adjust the distance between the Partial Reflectors to obtain a maximum meter reading. Record, d₁, the distance between the reflectors.



- **d**₁ = _____.
- 4. While watching the meter, slowly move one Reflector away from the other. Move the Reflector until the meter reading has passed through at least 10 minima and returned to a maximum. Record the number of minima that were traversed. Also record d_2 , the new distance between the Reflectors.
- Minima traversed = _____.
- **d**₂ = _____.
- 5. Use your data to calculate λ , the wavelength of the microwave radiation.
- λ = _____.
- 6. Repeat your measurements, beginning with a different distance between the Partial Reflectors.
- $\mathbf{d}_1 =$ _____. Minima traversed = _____.
- **d**₂ = _____. λ = _____.

Questions

- 1. What spacing between the two Partial Reflectors should cause a minimum signal to be delivered to the Receiver?
- 2. In an optical Fabry-Perot interferometer the interference pattern usually appears as a series of concentric rings. Do you expect such a pattern to occur here? Why? Check to see if there is one.

Experiment 9: Michelson Interferometer

Equipment Needed:

Item	Item
Transmitter	Goniometer
Receiver	Component Holders (2)
Partial Reflectors (1)	Reflectors (2)
Fixed Arm Assembly	Rotating Component Holder

Introduction

Like the Fabry-Perot interferometer, the Michelson interferometer splits a single wave, then brings the constituent waves back together so that they superpose, forming an interference pattern. Figure 9.1 shows the setup for the Michelson interferometer. A and **B** are Reflectors and **C** is a Partial Reflector. Microwaves travel from the Transmitter to the Receiver over two different paths. In one path, the wave passes directly through **C**, reflects back to **C** from **A**, and then is reflected from **C** into the Receiver. In the other path, the wave reflects from **C** into **B**, and then back through **C** into the Receiver.

If the two waves are in phase when they reach the Receiver, a maximum signal is detected. By moving one of the Reflectors, the path length of one of the waves changes, thereby changing its phase at the Receiver so a maxima is no longer detected. Since each wave passes twice between a Reflector and the Partial Reflector, moving a Reflector a distance $\lambda/2$ will cause a complete 360-degree change in the phase of one wave at the Receiver. This causes the meter reading to pass through a minimum and return to a maximum.

Purpose

In this experiment you will investigate the Michelson Interferometer.

Procedure

- 1. Arrange the equipment as shown in Figure 9.1. Plug in the equipment and adjust the Receiver for an easily readable signal.
- 2. Slide Reflector A along the Goniometer arm and observe the relative maxima and minima of the meter deflections.
- 3. Set Reflector A to a position which produces a maximum meter reading. Record, x_1 , the position of the Reflector on the Goniometer arm.
- **x**₁ = _____.
- 4. While watching the meter, slowly move Reflector A away from the Partial Reflector. Move the Reflector until the meter reading has passed through at least 10 minima and returned to a maximum. Record the number of minima that were traversed. Also record \mathbf{x}_2 , the new position of Reflector A on the Goniometer arm.
- Minima traversed = _____.
- $\mathbf{x}_2 =$ _____.
- 5. Use your data to calculate λ , the wavelength of the microwave radiation.
- λ = _____.
- 6. Repeat your measurements, beginning with a different position for Reflector A.







Experiment Guide

- Minima traversed = _____
- x₂ = _____.
- λ = _____.

Question

1. You have used the interferometer to measure the wavelength of the microwave radiation. If you already knew the wavelength, you could use the interferometer to measure the distance over which the Reflector moved. Why would an optical interferometer (an interferometer using visible light rather than microwaves) provide better resolution when measuring distance than a microwave interferometer?

An Idea for Further Investigation

Place a cardboard box between the Partial Reflector and Reflector A. Move one of the reflectors until the meter deflection is a maximum. Slowly fill the box with styrene pellets while observing the meter deflections. On the basis of these observations, adjust the position of Reflector A to restore the original maximum. Measure the distance over which you adjusted the reflector. Also measure the distance traversed by the beam through the pellets. From this data, can you determine the styrene pellets' index of refraction at microwave frequencies? The wavelength of electromagnetic radiation in a material is given by the relationship $\lambda = \lambda_0/\mathbf{n}$; where λ is the wavelength, λ_0 is the wavelength in a vacuum, and \mathbf{n} is the index of refraction of the material. Try boxes of various widths. You might also try filling them with a different material.

Experiment 10: Fiber Optics

Equipment Needed:

ltem	Item	
Transmitter	Goniometer	
Receiver	Styrene Pellets	
Tubular Plastic Bags		

Introduction

Light can propagate through empty space, but it can also propagate well through certain materials, such as glass. In fiber optics, a thin, flexible glass tube functions as a transmission line for light from a laser, much as a copper wire can function as a transmission line for electrical impulses. In the same way that variation of the electrical impulses can carry information through the copper wire (for example as a phone message), variation in the intensity of the laser light can carry information through the glass tube.

Purpose

In this experiment you will investigate the phenomenon of total internal reflection (as related to "fiber optics").

Procedure

- 1. Align the Transmitter and Receiver directly across from each other on the Goniometer, and adjust the Receiver controls for a readable signal.
- 2. Fill a tubular plastic bag with styrene pellets (tie the end or use a rubber band). Place one end of the bag in the Transmitter horn. What happens to the meter reading? Now place the other end in the Receiver horn. How does the intensity of the detected signal compare to the intensity when the bag is not used?
- **3.** Remove the plastic bag and turn the Rotatable Goniometer arm until no meter deflection appears. Place one end of the bag in the Transmitter horn, the other in the Receiver horn. Note the meter reading.
- 4. Vary the radius of curvature of the plastic bag. How does this effect the signal strength? Does the signal vary gradually or suddenly as the radial curvature of the plastic bag changes? Find the radius of curvature at which the signal begins to drop significantly.

Questions

- 1. Check your textbook for information on Total Internal Reflection. Based on the radial curvature when the signal begins to show attenuation as it passes through the plastic bag, determine the angle of total internal reflection for the styrene pellets. Can you use this value to determine the index of refraction of the styrene pellets?
- 2. Would you expect the plastic bag filled with styrene pellets to work the same with radiation at optical frequencies? Why?

Experiment 11: Brewster's Angle

Equipment Needed:

Item	Item	
Transmitter	Goniometer	
Receiver	Rotating Table	
Polvethylene Panel		

Introduction

When electromagnetic radiation passes from one media into another, some of the radiation usually reflects from the surface of the new medium. In this experiment, you will find that the magnitude of the reflected signal depends on the polarization of the radiation. In fact, at a certain angle of incidence—known as Brewster's Angle—there is an angle of polarization for which no radiation will be reflected. (Check your textbook for more information on Brewster's Angle.).

Purpose

In this experiment you will investigate the Brewster's Angle.

Procedure

- 1. Arrange the equipment as shown in Figure 11.1, setting both the Transmitter and the Receiver for horizontal polarization (90°).
- 2. Adjust the Panel so the angle of incidence of the microwave from the Transmitter is 20°. Rotate the Goniometer arm until the Receiver is positioned where it can detect the maximum signal reflected from the Panel. Adjust the Receiver controls for a mid-scale reading, and record the meter reading in Table 11.1.



Angle	Meter Reading (Horizontal Polarization)	Meter Reading (Vertical Polarization)
20°		
25°		
30°		
35°		
40°		
45°		
50°		
55°		
60°		
65°		
70°		
75°		

- 3. Without changing the angles between the transmitted beam, the Polyethylene Panel, and the Receiver, rotate both the Transmitter and the Receiver horns so they align for vertical polarization (0°) . Record the new meter reading in the table.
- 4. Repeat steps 2 and 3, setting the angle of incidence to each of the values shown in the table below. At each point set the Transmitter and Receiver for horizontal polarization and record the meter reading; then set them for vertical polarization and record that reading as well.
- 5. Plot a graph of "Meter Reading" versus "Angle of Incidence". Plot both the vertical and horizontal polarizations on the same graph. Label Brewster's Angle—the angle at which the horizontally polarized wave does not reflect.

Questions

- 1. Explain how Polaroid sun-glasses can be used to reduce the glare caused by the sun setting over a lake or the ocean. Should the glasses be designed to block vertically or horizontally polarized light?
- 2. Could you use the microwave apparatus to locate Brewster's Angle by examining the transmitted wave rather than the reflected wave? How?

Experiment 12: Bragg Diffraction

Equipment Needed:

Item	Item	
Transmitter	Goniometer	
Receiver	Rotating Table	
Cubic Lattice		

Introduction

Bragg's Law provides a powerful tool for investigating crystal structure by relating the interplanar spacings in the crystal to the scattering angles of incident x-rays. In this experiment, Bragg's Law is demonstrated on a macroscopic scale using a cubic "crystal" consisting of 10-mm metal spheres embedded in an ethafoam cube.

Before performing this experiment, you should understand the theory behind Bragg Diffraction. In particular, you should understand the two criteria that must be met for a wave to be diffracted from a crystal into a particular angle. Namely, there is a plane of atoms in the crystal oriented with respect to the incident wave, such that:

- 1. The angle of incidence equals the angle of reflection, and
- 2. Bragg's equation, $2d\sin\theta = n\lambda$, is satisfied; where d is the spacing between the diffracting planes, θ is the grazing angle of the incident wave, **n** is an integer, and λ is the wavelength of the radiation.

Purpose

In this experiment you will investigate Bragg Diffraction.

Procedure

1. Arrange the equipment as shown in Figure 12.1.



2. Notice the three families of planes indicated in Figure 12.2. (The designations (100), (110), and (210) are the Miller indices for these sets of planes.) Adjust the Transmitter and Receiver so that they directly face each other. Align the crystal so that the (100) planes are parallel to the incident microwave beam. Adjust the Receiver controls to provide a readable signal. Record the meter reading.

Experiment Guide

- **3.** Rotate the crystal (with the rotating table) one degree clockwise and the Rotatable Goniometer arm two degrees clockwise. Record the grazing angle of the incident beam and the meter reading. (The grazing angle is the complement of the angle of incidence. It is measured with respect to the plane under investigation, NOT the face of the cube; see Figure 12.3.)
- 4. Continue in this manner, rotating the Goniometer arm two degrees for every one degree rotation of the crystal. Record the angle and meter reading at each position. (If you need to adjust the INTEN-SITY setting on the Receiver, be sure to indicate that in your data.)



- 5. Graph the relative intensity of the diffracted signal as a function of the grazing angle of the incident beam. At what angles do definite peaks for the diffracted intensity occur?
- 6. UUse your data, the known wavelength of the microwave radiation (2.85 cm), and Bragg's Law to determine the spacing between the (100) planes of the Bragg Crystal. Measure the spacing between the planes directly, and compare with your experimental determination.
- 7. If you have time, repeat the experiment for the (110) and (210) families of planes.

Questions

- 1. What other families of planes might you expect to show diffraction in a cubic crystal? Would you expect the diffraction to be observable with this apparatus? Why?
- 2. Suppose you did not know beforehand the orientation of the "inter-atomic planes" in the crystal. How would this affect the complexity of the experiment? How would you go about locating the planes?

The Bragg Diffraction Experiment was developed by Dr. Harry Meiners of Rensselaer Polytechnic Institute.



Teacher's Guide

Experiment 1: Introduction to the Microwave Optics System

Notes on the Procedure

- 1. The meter reading does not vary with distance in an entirely predictable way, since the microwaves form standing waves between the transmitter and receiver at certain distances. In addition, the meter is not directly related to either the electric field or the intensity of the incident beam. The meter is useful for measuring relative intensity at a constant distance, polarization, and so on.
- **2.** The meter reading oscillates as the distance is decreased.
- 3. The presence of a reflector increases the meter reading.
- 4. The receiver detects no signal when the transmitter and receiver are at 90° to each other.



There is no significant difference between the output distributions in the horizontal and vertical orientations





Answers to Questions

- 1. The meter reading is not proportional to either the electric field or the intensity.
- 2. The meter reading is not proportional to either the electric field or the intensity.
- 3. The transmitter output is more plane wave than spherical wave, but it has characteristics of both.



Experiment 2: Reflection

Notes on the Procedure

5.

Table 2.1: Data	
Angle of Incidence	Angle of Reflection
20°	23°
30°	31°
40°	41°
50°	52°
60°	63°
70°	85°
80°	78°
90°	70°

The last three points are suspect due to the spread of the output pattern of the transmitter. See Experiment 1, part 8.

Answers to Questions

- 1. The angle of incidence equals the angle of reflection. This does hold for all angles, although it is not clear in this experiment due to the spread in the output pattern.
- 2. Some of the wave appeared to reflect into different angles; particularly when the angle of incidence was 70° or 90°. This is actually a diffraction effect, not reflection.
- 3. The transmitter does not produce a perfect plane wave, and this does affect the results.

Answers to Questions for Additional Experimentation

- 1. Intensity of the reflection does vary with the angle of incidence; from this we can deduce that the reflector is not 100% efficient.
- 2. In general, conductors will reflect the microwaves much better than non-conductors.

3.

Experiment 3: Standing Waves - Measuring Wavelengths

Notes on the Procedures

Method A

Antinodes	Distance (cm)	Wavelength (cm)
10	13.3	2.66
15	20.5	2.73

Average wavelength: 2.70 cm

Frequency: 1.11 x 10¹⁰ Hz

(Fewer points were taken due to the limited resolution of this method.)



Method B

Antinodes	Distance (cm)	Wavelength (cm)
5	7.1	2.84
10	14.1	2.82
15	20.0	2.67
20	27.5	2.89

Average wavelength: 2.81 cm

Frequency: 1.07 x 10¹⁰ Hz

Answers to Questions

The value obtained by the first method was 5.6% off, and the second was 1.5% off. If it is possible to take more data points on the second method, you may get better results.

Experiment 4: Refraction Through a Prism

Notes on the Procedure

- 1. The empty foam prism absorbs the radiation by a very slight amount.
- 3. $\theta = 7^{\circ} (\pm 1^{\circ})$
- $\mathbf{4.} \quad \mathbf{\theta}_1 = 22^{\circ}$
 - $\theta_2 = 29^\circ$
- 5. Our experimental value was: $\mathbf{n}_1 = 1.3 \pm 0.05$
- 6. Our experimental value was: $\mathbf{n}_1 = 1.3 \pm 0.05$

Answers to Questions

- 1. This assumption is valid. According to Snell's law, if the angle of incidence is zero, the angle of refraction is zero also.
- 2. The index for a solid styrene prism would be higher, due to the greater "optical" density of the solid material.

General Notes

The prism mold may be filled with other materials as well. We used water for one such test. The water absorbs most of the microwave energy (this is how a microwave oven works) but enough gets through that it may be measured on the most sensitive scale of the receiver. We found that $\mathbf{n} = 1.4 \pm 0.05$.

Experiment 5: Polarization

Notes on the Procedure

2.





NOTE: There is a consistent "glitch" in the data at a polarization angle of about 40 and 140 degrees which is not entirely explained by the non-linearity of the receiver. (This glitch is also present when the polarizer slits are used in part 5 of this lab.)

- 3. Continued rotation of the receiver results in duplications of the pattern above.
- 5.



6. The meter reading is zero when the polarizer slits are oriented vertically or horizontally. When the slits are at 45°, the meter reads about 30% of its maximum value for that distance.

Answers to Questions

- 1. The meter reading more closely matches the intensity than it does the electric field.
- 2. The meter reading more closely matches the intensity than it does the electric field.



3. The polarizer transmits only the component of the wave parallel to the polarizer.



4. When the transmitter and polarizer are at 90°, the wave is completely blocked. Placing a polarizer at 45°, however, introduces a component of the wave parallel to the receiver so that some of the wave is then picked up.



Experiment 6: Double-Slit Interference

Answers to Questions



- 1. For the 7.6 cm spacing, the maxima should occur at 22° and 48°.
- 2. For the 10.6 cm spacing, the maxima should be at 15° , 33° , and 54° .
- These theoretical values are closely matched by the experimental data. The theory assumes that the distance between slit plates and receiver are large compared to the slit spacing and wavelength. This requirement is barely met in the experimental setup used, and could cause trouble in some situations.
- 3. The single-slit pattern (see experiment 7) acts as an upper limit to the multiple-slit pattern from this experiment.

General Notes

1. The position of the experimenter has a definite effect on the measurements in this experiment. Experiment to find just how much effect there is with your particular setup, and then take your data accordingly.

2. Single-slit diffraction may also be attempted on this apparatus, though we don't recommend it. The distances are too short, relative to the wavelength; so the analysis requires the Fresnel/Kirchoff approach instead of the Fraunhoffer approximation. Even with the Fresnel approach, the "fringes" are too small to be seen adequately.



• Basically, it doesn't work well at all.

Experiment 7: Lloyd's Mirror

Notes on the Procedure

Fringe	h (mm)	Path length (mm)	Wavelength (mm)
1	85	720	
2	138	752	32.10
3	171	779	26.63
4	200	806	27.15
5	227	834	28.11
6	251	861	27.06
7	276	890	28.83
8	298	919	29.13
9	317	944	25.08
10	338	973	28.69

Average: 27.59 mm

Standard deviation: 1.37

An alternate method is to graph path length versus fringe number and take the slope of the graph. This slope will be the wavelength.

Answer to the Question

• It simplifies the calculations.





Experiment 8: Fabry-Perot Interferometer

Notes on the Procedure

1. through 4.

For best results, do not move the reflector closest to the transmitter. There are actually two standing wave patterns that may form: one between the transmitter and first reflector, and one between the two reflectors. (There may also be others, such as between the second reflector and the receiver or the second reflector and the transmitter; but these will be negligible.) Moving the first reflector will change the amplitude of the wave coming into the region between the reflectors, and thus give erroneous results.

n	Second plate	Distance	D	n	Second plate	Distance	D
1	54.9	20.3		12	39.4	35.8	1.4
2	53.6	21.6	1.3	13	38.0	37.2	1.4
3	52.2	23.0	1.4	14	36.6	38.6	1.4
4	50.9	24.3	1.3	15	35.1	40.1	1.5
5	49.4	25.8	1.5	16	33.8	41.4	1.3
6	47.9	27.3	1.5	17	32.3	42.9	1.5
7	46.5	28.7	1.4	18	30.9	44.3	1.4
8	45.0	30.2	1.5	19	29.4	45.8	1.5
9	43.6	31.6	1.4	20	27.9	47.3	1.5
10	42.1	33.1	1.5	21	26.5	48.7	1.4
11	40.8	34.4	1.3	22	25.0	50.2	1.5
	•					Average	1.42

First plate: 75.2 mm

5. and 6.

2.85 cm

NOTE: An alternate method of analysis is to make a graph of distance versus fringe number and take the slope of the line to find the wavelength.



Answers to Questions

- 1. Minima will occur when the spacing is $n\lambda/4$, where n is an odd integer.
- 2. Normally, just such a pattern would occur; in this case, however, the reflectors are too small in relation to the wavelength used; so the next "ring" is located beyond the edge of the reflectors and may not be seen.

Experiment 9: Michelson Interferometer

Notes on the Procedure

1. through 4.

Best results are obtained when the mirrors are both a significant distance from the central beam-splitter. If either mirror is too close to the center, the maxima splits into two peaks due to secondary interference effects.

n	Refl. Pos.		n	Refl. Pos.	
1	75.0		12	90.9	1.5
2	76.4	1.4	13	92.3	1.4
3	77.9	1.5	14	93.8	1.5
4	79.4	1.5	15	95.2	1.4
5	80.8	1.4	16	96.6	1.4
6	82.2	1.4	17	98.1	1.5
7	83.7	1.5	18	99.5	1.4
8	85.1	1.4	19	101.0	1.5
9	86.6	1.5	20	102.4	1.4
10	88.0	1.4	21	104.0	1.6
11	89.4	1.4		Average	1.45

From this it can be calculated that the wavelength is 2.80 cm

NOTE: An alternate method of analysis is to plot the reflector position versus fringe number. The slope of this line will be half the wavelength.

By this method, the wavelength is calculated as 2.89 cm.

Answer to Question

1. The limit of resolution for distance measurements with a Michelson interferometer is roughly 1/4 the wavelength of the light used. Thus with these microwaves, distance changes of about 7 mm can be measured. With a visible-light interferometer and a wavelength of 633 nm (HeNe laser light), distance changes of only 158 nm can be measured.

Experiment 10: Fiber Optics

Notes on the Procedure

- 2. The meter reading with the "optical fiber" bag in place can more than double the unobstructed meter reading. This is because the bag prevents the normal spreading of the beam, and directs all the microwave radiation into the receiver horn.
- 3. The bag can direct the beam through a full 90° turn without measurable attenuation, if the curvature is gradual.
- 4. The signal begins to be attenuated with a radius of curvature of about 5 cm, and drops off rather suddenly from there. It is difficult to get consistent results, though.

Answers to Questions

1. Theoretically one can use the radius at which the microwaves begin to "leak" to determine the index of refraction of the material. In reality, this is quite difficult. Our values for n using this method range from 1.1 to 1.4.



2. The styrene pellets are very small, compared to the microwave wavelength. Compared to visible wavelengths, the pellets are enormous. Because of this size difference, optical radiation is scattered by the pellets and microwave radiation is transmitted.

Additional Idea

This apparatus may also be used as a demonstration of how the plane of polarization can be rotated by multiple reflections. Rotate the transmitter or receiver 90° to each other, so that no signal gets through. Now put the bag of pellets between the two. If the bag is held straight, there will be a zero meter reading; but if the bag is curved into a spiral, there will be a non-zero reading.

Experiment 11: Brewster's Angle

Notes on the Procedure

1. through 5.



Answers to Questions

- 1. Glare off water from a low source is primarily horizontal in polarization, so sunglasses should be designed to block horizontally polarized light.
- 2. Theoretically, one could do this by finding the point at which there was no transmission of vertically polarized light. We have not been able to get good results in this experiment, however.

General Notes

- 1. The index of refraction of the polyethylene is 1.5. (Calculated from the dielectric constant at 10 GHz)
- 2. One must be careful on this experiment to note that there are actually two effects which are being measured. In addition to the Brewster's angle reflection, there is a certain amount of interference between the front- and rear-surface reflections from the polyethylene. This interference causes extra peaks at 13.8° and 67.2°, and a local minimum at 43.6°. The interference is notable enough that you may want to demonstrate the interference instead of the Brewster's angle.

Experiment 12: Bragg Diffraction

Notes on the Procedure

1. The polarization angle of the transmitter and receiver (horizontal or vertical) does not matter.

5.



Peaks occur at 18° , 24° , and 45° . These correspond to plane spacings of 4.6cm, 3.4cm, and 4.0 cm (n = 2 for the 45° peak) The actual spacing is 3.8 cm The first peak is apparently a reflection off a different plane than the one we're measuring.

6. 110 Plane:



The peak at 29° gives a plane spacing of 2.9 cm; the actual spacing is 2.7 cm

Answers to Questions

- 1. Other families of planes would be the 111 plane, the 101 plane, and so on. These would be difficult to observe with this apparatus due to the small size of our "crystal".
- 2. Not knowing the orientation of the interatomic planes of the crystal would increase the complexity of the analysis. We could orient the crystal so that there was maximum transmission; this would indicate to us a 100 plane. From there, we could try different angles until we had enough data to assemble a likely picture of the atomic spacing in the crystal.
- NOTE: Error of 10 to 15% is reasonable for this experiment.

