
#### Instruction Sheet for the PASCO Model SE-7256

## Motion Sensor Guard



Introduction

The SE-7256 Motion Sensor Guard is designed to protect the PASCO Motion Sensor from falling objects. When placed over the motion sensor, the wires of the guard are within the 15 cm "dead zone" of the sensor. Therefore, the wires do not interfere with data collection for position, velocity or acceleration measurements.


#### **Setup Instructions**

- 1. Set the Motion Sensor on a flat surface
- 2. Rotate the sensor head to point upwards
- 3. Place the Motion Sensor Guard over the Motion Sensor as shown.
- 4. Be sure to position the guard so the intersections of the wires are not directly above the sensing element. This will minimize any interference by the guard.

**Note:** The Motion Sensor Guard is designed to protect the sensor from light objects such as a foam ball, coffee filters, or a tennis ball. Heavy objects may break the guard and damage the motion sensor.



Position the guard so that the intersections of the wires are not directly above the sensing element.







Motion Sensor Guard 012-08901B

### Warranty

For a description of the product warranty, see the warranty description in the PASCO catalog or on the PASCO web site at www.pasco.com.

# Copyright

The PASCO scientific Model SE-7256 Motion Sensor Guard instruction sheetl is copyrighted and all rights reserved. However, permission is granted to non-profit educational institutions for reproduction of any part of the instruction sheet providing the reproductions are used only for their laboratories and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

## **Technical Support**

Address: PASCO scientific

10101 Foothills Blvd. Roseville, CA 95747-7011

Phone: (916) 786-3800 or toll-free in the US

(800)772-8700

FAX: (916)786-7565

email: techsupp@PASCO.com

web: www.pasco.com

