


Micro Stir Bar

PS-2565



| Included Equipment | Part Number |
|-------------------------------|------------------------------------|
| Micro Stir Bar (5 pieces) | PS-2565 |
| Additional Equipment Required | |
| Stir Plate | SE-7700 or similar |
| Sensor with probe | See PASCO catalog or www.pasco.com |

Introduction

The PASCO Micro Stir Bar maintains a constant flow of solution over the end of a probe. PASCO part PS-2565 is package of five stir bars.

The stir bar attaches to the end of a probe. Use a standard stir plate (such as PASCO model SE-7700) to drive the stir bar.

The stir bar is compatible with cylindrical probes of about 13 mm in diameter (such as *ScienceWorkshop* and PASPORT pH, conductivity, and dissolved oxygen probes).

Set-up

1. Place the stir bar over the end of the probe and press it into place.



- 2. Place a beaker on a stir plate.
- **3.** Put the probe into the beaker. Use a stand and clamp to hold the probe.
- **4.** Turn on the stir plate to make the stir bar spin.

Technical Support

For assistance with any PASCO product, contact PASCO at:

Address: PASCO scientific

10101 Foothills Blvd. Roseville, CA 95747-7100

Phone: 916-786-3800 (worldwide)

800-772-8700 (U.S.)

Fax: (916) 786-7565 Web: www.pasco.com Email: support@pasco.com

Limited Warranty For a description of the product warranty, see the PASCO catalog.

Copyright The PASCO scientific 012-10251A *Micro Stir Bar Instruction Sheet* is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited. **Trademarks** PASCO, PASCO scientific, PASPORT, and

ScienceWorkshop are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit www.pasco.com/legal.

800-772-8700 www.pasco.com