


Super-Flex I-Beam

ME-8987



Included Equipment	Part Number
Super-Flex I-Beam	ME-8987

Introduction

The model ME-8987 Super-Flex I-Beam illustrates beam stiffness and deflection. Because it flexes more than a metal or wooden beam, students can easily observe that its vertical stiffness is greater than its horizontal stiffness. Use the Super-Flex I-Beam to qualitatively demonstrate how a beam bends under various load and support configurations.

Grids on the top and one side of the I-beam show the amount of deflection in a qualitative way. The grid on top is approximately 25.4 cm (10 in) long and 3.8 cm (1.5 in) wide with 1.2 cm by 1.2 cm (0.5 in by 0.5 in) squares. The grid on the side is approximately 12 cm long and 2 cm wide with 1 cm by 1 cm squares.

Demonstrations

Bending

Hold the I-beam in your hands and bend it in the “upright” and “side” directions as illustrated. Hold the I-beam with the top grid facing the students and demonstrate the distortion as the beam is bent (see Figure 2a). In which direction is it easier to bend? For the same effort, in which direction does it bend more?




Figure 1: Illustration of I-beam bent in upright and side directions

Buckling

Bend the I-beam in the upright position and continue to bend it until the beam buckles (see Figure 2b). Another way to show this mode of failure is to stand the I-beam on its end and press down on the top end until the column buckles.

Torsion

Twist the I-beam to demonstrate its lack of strength under torsion (see Figure 2c).




Figure 2

Technical Support

For assistance with any PASCO product, contact PASCO at:

Address: PASCO scientific
10101 Foothills Blvd.
Roseville, CA 95747-7100

Phone: 916-786-3800 (worldwide)
800-772-8700 (U.S.)

Fax: (916) 786-7565

Web: www.pasco.com

Email: support@pasco.com

Limited Warranty

For a description of the product warranty, see the PASCO catalog.

Copyright

The PASCO scientific 012-10603A *Super-Flex I-Beam Instruction Sheet* is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

Trademarks

PASCO, PASCO scientific, and ScienceWorkshop are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit www.pasco.com/legal.