


Wireless Centripetal Force Accessory

ME-8094



1	Centripetal Force Accessory Platform	4	Sensor Mounting Rod
2	"Tee" Tab	5	Thumbscrew (X2)
3	Square Nut (1/4-20 by 7/32 inch)	6	Bushing (X2)

Hardware


Included Equipment	Part Number		
Wireless Centripetal Force Accessory	ME-8094		
Recommended Equipment	Part Number		
Wireless Force-Acceleration Sensor	PS-3202		
Rotating Platform	ME-8951		
Hanger and Mass Set	ME-8979		
Physics String (SE-8050) or Thread (ME-9875)			

Introduction

The Wireless Centripetal Force Accessory is a low-friction sliding mass holder that connects to a Wireless Force Acceleration Sensor. When installed on a Rotating Platform it supports the direct measurement of centripetal force and acceleration. It can be used with the Wireless Force Acceleration Sensor to explore the relationships among the variables involved in centripetal acceleration and centripetal force.

Vary the mass using the holed masses in the Mass and Hanger Set and vary the radius by simply adjusting the string length. Vary the angular velocity by spinning the Rotating Platform faster or slower. The Wireless Force Acceleration Sensor measures force along the axis between the metal plate in the front and the thumbscrew at the end. It measures acceleration in three directions (X-Y-Z). When mounted on the Rotating Platform, the sensor measures angular velocity around the Z-axis.

Setup



Slide the "Tee Tab" on the bottom of the Wireless Centripetal Force Accessory into the "Tee" slot on the top of the Rotating Platform. Loosen the square nut on the Sensor Mounting Rod and slide it into the same "Tee" slot. Put the Hook on the end of the Wireless Force Acceleration Sensor and mount the sensor on the Sensor Mounting Rod. Position the sensor so that the end of the hook is at the zero mark on the top of the Rotating Platform. Tighten the Sensor Mounting Rod in position.


Attach one end of a 30 cm piece of string or thread to the hook on the sensor. Attach the other end to a thumbscrew on the Wireless Centripetal Force Accessory platform (See Details).

Details

Loosen the thumbscrew. Put the string or thread under the bushing as shown. Tighten the thumbscrew. To change the length of the string or thread, loosen the thumbscrew.

Usage

See the PASCO web site at www.pasco.com. Go to the web page for the Wireless Centripetal Force Accessory and click "Video" to see a short demonstration.



Technical Support

For assistance with any PASCO product, contact PASCO at:

Address: PASCO	scientific
10101 F	Foothills Blvd.
Rosevil	le, CA 95747-7100

- Phone: 916-786-3800 (worldwide) 800-772-8700 (U.S.)
- Fax: (916) 786-7565
- Web: www.pasco.com
- Email: support@pasco.com

More Information

For more information about the latest revision of this Instruction Sheet, visit:

www.pasco.com/manuals

and enter the Product Number.

For information about the Wireless Centripetal Force Accessory or any PASCO product, what software to use, and what other accessories are available, check the PASCO web site:

www.pasco.com

Warranty, Copyright, and Trademarks

Limited Warranty For a description of the product warranty, see the PASCO catalog. Copyright The PASCO scientific 012-14793A *Wireless Centripetal Force Accessory Instruction Sheet* is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited. Trademarks PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit www.pasco.com/legal.