


# Wireless Heart Rate Sensor Hand Grips | PS-3565

**PASCO**

Quick Start Guide  
012-16973A



## Required Equipment


The Wireless Heart Rate Sensor Hand Grips must be used with the Wireless Heart Rate Sensor Module (PS-3566). The module and Wireless Heart Rate Sensor Hand Grips are both included with the PS-3206.

## For Educational Use Only

**Note:** This is not a medical device. It is designed for educational use only and should not be used in any medical process such as life support or patient diagnosis. It is also not intended for use in graduate research or industry including industrial control or any type of industrial testing.

## Using the Hand Grips

1. Remove the protective tape from the metallic strips on the front and back of each grip.
2. Snap the Heart Rate Module onto the hand grip connectors.



3. Grasp the hand grips firmly (but not too tightly) in the palms of each hand so that one of the metallic strips is in the center of the palm.

## Care and Maintenance

Remove the Wireless Heart Rate Monitor from the Hand Grips when the sensor is not in use. Clean the hand grips after each use with a damp cloth or paper towel. The hand grips can be cleaned with anti-bacterial cleaners. Do not submerge the hand grips or the Wireless Heart Rate Module in water.

## Troubleshooting

Erratic measurements or large reading fluctuations usually indicate signal loss. To avoid signal loss:

- Check the position of the hands on the hand grips. Make sure both palms are touching the metal strips of each hand grip.
- Avoid close proximity to computer monitors, magnets, power cords, television sets, or other appliances that emit electromagnetic radiation.

If necessary, restart the software.

## Experiments

Print-ready experiment worksheets are available to download from the PASCO website. Go to [pasco.com/freelabs](http://pasco.com/freelabs) and enter **PS-3206** in the **Part No.** field to access all experiments related to the Wireless Heart Rate Sensor Hand Grips.

## Technical Support

Contact PASCO Technical Support for assistance with this or any other PASCO product.

Address	PASCO scientific 10101 Foothills Blvd. Roseville, CA 95747-7100
Phone	1-800-772-8700 (USA) 1-916-462-8384 (International)
Chat	<a href="http://www.pasco.com">www.pasco.com</a>
email	<a href="mailto:support@pasco.com">support@pasco.com</a>

## Regulatory Information

### Warranty, Copyright, and Trademarks

#### Limited Warranty

For a description of the product warranty, see the Warranty and Returns page at [www.pasco.com/legal](http://www.pasco.com/legal).

#### Copyright

This document is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

#### Trademarks

PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit [www.pasco.com/legal](http://www.pasco.com/legal).

### Product End of Life Disposal Instructions



This electronic product is subject to disposal and recycling regulations that vary by country and region. It is your responsibility to recycle your electronic equipment per your local environmental laws and regulations to ensure that it will be recycled in a manner that protects human health and the environment. To find out where you can drop off your waste equipment for recycling, please contact your local waste recycle or disposal service, or the place where you purchased the product.

The European Union WEEE (Waste Electronic and Electrical Equipment) symbol and on the product or its packaging indicates that this product must not be disposed of in a standard waste container.

### FCC Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

### CE Statement

This device has been tested and found to comply with the essential requirements and other relevant provisions of the applicable EU Directives.