


Cart Masses

ME-6757A


Introduction

The Cart Masses are a pair of 250 g (0.25 kg) masses designed to be used with any PASCO dynamics cart. These masses provide an easy way to temporarily increase the mass of a cart during an experiment.



Usage

The Cart Masses are designed to be mounted in the tray of any PASCO dynamics cart, such as the Smart Cart (ME-1240 or ME-1241), the PAScar (ME-6933 or ME-6934), the Plunger Cart (ME-9430), or the Collision Cart (ME-9454). Each mass has approximately the same mass as a single plastic cart (including the Smart Cart) and half the mass of a metal cart. The cart tray can accommodate two masses. Additional masses can be added by stacking the bars, as shown below. The built-in grooves on the underside of the bars keep the masses in place.



Technical support

Need more help? Our knowledgeable and friendly Technical Support staff is ready to answer your questions or walk you through any issues.

- Chat [pasco.com](https://www.pasco.com)
- Phone 1-800-772-8700 x1004 (USA)
+1 916 462 8384 (outside USA)
- Email support@pasco.com

Limited warranty

For a description of the product warranty, see the Warranty and Returns page at www.pasco.com/legal.

Copyright

This document is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

Trademarks

PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit www.pasco.com/legal.