


# Economy CASTLE Kit



The Economy CASTLE Kit includes four complete CASTLE Kits. A complete kit is shown in the figure. The storage box is *not* included with the Economy CASTLE Kit.

## 1 Light Bulb Sockets

- 16 included (4 per CASTLE Kit)
- Replacement part: EM-8630 (10 pack)

## 2 Battery Holder

- 4 included (1 per CASTLE Kit)
- Accepts four D-cell batteries (not included)

## 3 Compass

- 4 included (1 per CASTLE Kit)
- Replacement part: EM-8631A (5 pack)

## 4 Capacitor, 0.025 F

- 4 included (1 per CASTLE Kit)
- Replacement part: EM-8632 (2 pack)

## 5 Alligator Clip Leads

- 40 included (10 per CASTLE Kit)
- Replacement part: EM-8634 (Set of 10)

## 6 Light bulb, #48, 2 V, 60 mA, oblong

- 25 included (6 per CASTLE Kit)

## 7 Light bulb, #14, 2.5 V, 300 mA, round

- 25 included (4 per CASTLE Kit)

## 8 Resistor, 10 $\Omega$ , 1 W (not pictured)

- 16 included (4 per CASTLE Kit)

## Experiments

Copy-ready experiment worksheets for classroom use are available on the PASCO website. Download the experiments for free at [www.pasco.com/resources/lab-experiments](http://www.pasco.com/resources/lab-experiments).

## Technical Support

Contact PASCO Technical Support for assistance with this product or any other PASCO product.

Address	PASCO scientific 10101 Foothills Blvd. Roseville, CA 95747-7100
Phone	1-800-772-8700 (USA) 1-916-462-8384 (International)
Chat	<a href="http://www.pasco.com">www.pasco.com</a>
email	<a href="mailto:support@pasco.com">support@pasco.com</a>

## CE Statement

This device has been tested and found to comply with the essential requirements and other relevant provisions of the applicable EU Directives.

## Product End of Life Disposal Instructions



This electronic product is subject to disposal and recycling regulations that vary by country and region. It is your responsibility to recycle your electronic equipment per your local environmental laws and regulations to ensure that it will be recycled in a manner that protects human health and the environment. To find out where you can drop off your waste equipment for recycling, please contact your local waste recycle or disposal service, or the place where you purchased the product.

The European Union WEEE (Waste Electronic and Electrical Equipment) symbol and on the product or its packaging indicates that this product must not be disposed of in a standard waste container.

## Warranty, Copyright, and Trademarks

**Limited Warranty:** For a description of the product warranty, see the Warranty and Returns page at [www.pasco.com/legal](http://www.pasco.com/legal).

**Copyright:** This document is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

**Trademarks:** PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit [www.pasco.com/legal](http://www.pasco.com/legal).