


# Large Rod Base

ME-8735

## Introduction

The Large Rod Base can be used whenever you need a sturdy base for a vertical support rod. Rods with an outside diameter of  $\frac{1}{4}$ " to  $\frac{1}{2}$ " (6.5 mm to 13 mm) may be placed in either the central mounting hole or the alternate mounting hole at the peak of the stand. The rods can be secured in place using either of the two clamp screws provided.

## Equipment included



- ❶ Large Rod Base (ME-8735)
- ❷ 2× leveling feet
- ❸ 2× clamp screws

## Usage tips

- The mass of the Large Rod Base is approximately 4 kg. If dropped, the base may be damaged, cause damage to other equipment, or cause injury. To prevent injury or damage, always carry the Large Rod Base with a firm grip and, whenever possible, set up the stand away from the edge of any table or work bench.
- To prevent the Large Rod Base from tipping over during experiments, be careful as you set up experiments. Always make sure the Large Rod Base can rest solidly in any experiment involving motion.

## Specifications and accessories

Visit the product page at [pasco.com/product/ME-8735](https://www.pasco.com/product/ME-8735) to view the specifications and explore accessories. You can also download support documents from the product page.

## Technical support

Need more help? Our knowledgeable and friendly Technical Support staff is ready to answer your questions or walk you through any issues.

- Chat [pasco.com](https://www.pasco.com)
- Phone 1-800-772-8700 x1004 (USA)  
+1 916 462 8384 (outside USA)
- Email [support@pasco.com](mailto:support@pasco.com)

## Limited warranty

For a description of the product warranty, see the Warranty and Returns page at [www.pasco.com/legal](https://www.pasco.com/legal).

## Copyright

This document is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

## Trademarks

PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit [www.pasco.com/legal](https://www.pasco.com/legal).