Egg Drop – **Master Materials and Equipment**

Equipment by Activity

NOTE: Italicized entries indicate items not available from PASCO. The quantity indicated is per student or group.

Act	Title	Materials and Equipment	Qty
1	Reading Graphs	Data Collection System	1
	Use a motion sensor to generate graph data for practice in interpreting graphs.	PASPORT Motion Sensor	1
		Meter stick	1
		Pencils	2
2	Acceleration and Gravity	Data Collection System	1
	Use a motion sensor to analyze how gravity affects the motion of falling objects of different size and mass.	PASPORT Motion Sensor	1
		Large base and support rod	1
		Right angle clamp	1
		Small rod	1
		Balance, 1-g resolution	1 per class
		No-bounce pad	1
		Three objects of different size and similar mass	1
		Three objects of different mass and similar size	1
		Meter stick	1
3	Force	Data Collection System	1
	Use a motion sensor and force sensor to analyze the relationships between force, mass, and acceleration.	PASPORT Motion Sensor	1
		PASPORT Force Sensor with hook	1
		Compact cart mass, 250-g	2
		Dynamics cart	1
		Dynamics track	1
		Super pulley with clamp	1
		Adjustable end stop	1
		Hanging mass, 40-g	1
		Balance, 1-g resolution	1 per class
		String	85 cm
4	Air Drag Use a motion sensor to analyze how the velocity of an object in free fall is related to the net force on the object and that air drag produces a force opposite the direction of fall.	Data Collection System	1
		PASPORT Motion Sensor	1
		Large base and support rod	1
		Small rod	1
		Right angle clamp	1
		Balance, 0.01-g resolution	1 per class
		Coffee filters	4
		Meter stick	1
5	Impact Force	Data Collection System	1
	Use a force sensor to analyze force during a collision and determine the best material to cushion and protect fragile objects during a collision.	PASPORT Force Sensor with rubber bumper	1
		Discover Collision Bracket with spring bumpers	1
		Force sensor cup	1
		Large base and support rod	4
		Dynamics cart	1
		Dynamics track with feet	1

Act	Title	Materials and Equipment	Qty
		Track rod clamp	1
		Different bumper materials	4
		Tape	1 roll