

# Welcome to the



With the addition of the Xplorer GLX your students will begin a new chapter in Hands-on, Minds-on science education.

We know you are anxious to start using your GLX, and the enclosed Quick Start Poster will get you up and running, even doing your first experiment. But before you start, just take a minute to read this page, it's the most important investment in your GLX.

## Your GLX will get even better. Keep it updated!

Although your GLX is the most powerful and versatile handheld science computer, it's only going to get better and better. New firmware features and enhancements are currently in development.

### Updating is easy:

- go to [www.pasco.com/glx](http://www.pasco.com/glx),
- click on the "Update" button and follow directions.
- Or, go to Chapter 5, "Firmware upgrade" in your Xplorer GLX manual.

Remember, GLX firmware updates are FREE, and with each update you get even more value from your GLX.

## Help us help you. Keep updated.

To keep you fully advised about:

- FREE enhancements to your GLX
- New Sensors, and
- New curricular materials

please let us notify you by e-mail. We promise not to visit your mailbox more than once a month, only give you worthwhile information, not to give your e-mail address to anyone, and to take you off the list if you should desire.

To provide us with your e-mail, please use the registration card enclosed, or go to [www.pasco.com/glx](http://www.pasco.com/glx) and register. **Do it now, it will ensure you get the most from your GLX.**

## A growing list of Compatible Printers.

In case you didn't know, you can print **Directly** from your GLX to selected H. P. printers. Yes, just plug the printer into the USB port and print graphs, tables, notes, etc, and with even greater resolution than the GLX screen. Below is a list of compatible printers. As new printers are added, they will be listed on the GLX web page. Before you buy a new printer for your GLX, check the latest list at: [www.pasco.com/glx](http://www.pasco.com/glx).

DeskJet 3820	DeskJet 5650
DeskJet 920c	DeskJet 940c
DeskJet 960c	DeskJet 995ck
DeskJet 6122	DeskJet 6840
PhotoSmart 7760	PhotoSmart 7660
PhotoSmart 7260	DeskJet 9300
DeskJet 6127	LaserJet 2300
DeskJet 5150	DeskJet 5740
DeskJet 840c	DeskJet 932c
DeskJet 950c	DeskJet 990cxi
PhotoSmart 7450	DeskJet 5850
DeskJet 6540	DeskJet 450
DeskJet 9650	Business InkJet 1100
Color InkJet cp 1700	PhotoSmart 7960
PhotoSmart 8450	

The PASCO logo consists of the word "PASCO" in a bold, italicized, sans-serif font. The letters are closely spaced and have a slightly slanted appearance.

10101 Foothills Blvd  
Roseville, CA 95747-7100  
1-800-772-8700

012-09307B