

PS-3587 WI Sensor Storage Tray for PS-3215 Colorimeter & Turbidity Sensors

Included Items	Part Number
Gratnells Tray Label	646-15891
Gratnells Tray Lid	699-15383
Gratnells Tray	699-271
Bottom Foam	740-15777
Vacuum Formed Tray	740-15812
Upper Foam Assembly	740-16072

Product Description

The PS-3587 is a custom storage system designed for the PS-3215 Wireless Colorimeter and Turbidity Sensor. The storage tray holds up to 10 sensors and accessories.

Step 1

Installing Label

Center the label side to side and place label about an $\frac{1}{4}$ " up from the bottom edge on the Gratnells Tray as shown in picture.

Step 2


Lower Sensor Group

Place the bottom foam piece into bottom of tray.

Place up to 5 sensors and 5 Cuvette Racks into the foam as shown.









800-772-8700 www.pasco.com

Step 3

Upper Sensor Group

Place the upper foam assembly on top of the lower sensors in tray.

Place up to 5 sensors and 5 Cuvette Racks into the upper foam assembly.

Step 4

Miscellaneous Item Tray

Next place the vacuum formed tray on top of the upper sensor tray. This tray is for miscellaneous items like the USB charging cables and the 100 NTU calibration bottles.

Step 5

Tray Lid

Finally, orient tray lid and snap lid onto tray as shown in picture.

Technical Support

For assistance with PASCO products, contact PASCO at:

Address: PASCO scientific

10101 Foothills Blvd.


Roseville, CA 95747-7100

Phone: 916-462-8384
Web: www.pasco.com
Email: support@pasco.com

Warranty, Copyright, and Trademarks

Limited Warranty For a description of the product warranty, see the PASCO catalog.

Copyright The PASCO scientific 012-0000A Item Name Instruction Sheet is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.





Trademarks PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit www.pasco.com/legal.

